

Fibroblast growth factor receptor 1

Catalog No: #AP83898



Package Size: #AP83898-1 50ug #AP83898-2 100ug #AP83898-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	FGFR-1,Basic fibroblast growth factor receptor 1,Tyrosine kinase receptor CEK1,CEK1
Accession No.	Uniprot ID:P21804
Uniprot	P21804
GeneID	396516;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FGFR1

Note: This product is for in vitro research use only