

Heat shock 70 kDa protein 1A

Catalog No: #AP83901



Package Size: #AP83901-1 50ug #AP83901-2 100ug #AP83901-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock 70 kDa protein 1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Heat shock 70 kDa protein 3,HSP70.3,Hsp68,Hsp70-3,Hsp70A1
Accession No.	Uniprot ID:Q61696
Uniprot	Q61696
GeneID	193740;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspa1a

Note: This product is for in vitro research use only