

# HMG box-containing protein 1

Catalog No: #AP83905



Package Size: #AP83905-1 50ug #AP83905-2 100ug #AP83905-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HMG box-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	HMG box transcription factor 1,High mobility group box transcription factor 1
Accession No.	Uniprot ID:Q2KJ34
Uniprot	Q2KJ34
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HBP1

Note: This product is for in vitro research use only