

## E3 ubiquitin-protein ligase HACE1

Catalog No: #AP83906



Package Size: #AP83906-1 50ug #AP83906-2 100ug #AP83906-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	E3 ubiquitin-protein ligase HACE1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	HECT domain and ankyrin repeat-containing E3 ubiquitin-protein ligase 1,HECT-type E3 ubiquitin transferase HACE1,RCJMB04_37k12
Accession No.	Uniprot ID:E1C656
Uniprot	E1C656
GeneID	421788;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HACE1

Note: This product is for in vitro research use only