

Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-3

Catalog No: #AP83911

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP83911-1 50ug #AP83911-2 100ug #AP83911-3 1mg

Description

Product Name	Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Transducin beta chain 3
Accession No.	Uniprot ID:Q61011
Uniprot	Q61011
GeneID	14695;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Gnb3

Note: This product is for in vitro research use only