

HMG box-containing protein 1

Catalog No: #AP83922



Package Size: #AP83922-1 50ug #AP83922-2 100ug #AP83922-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HMG box-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HMG box transcription factor 1,High mobility group box transcription factor 1
Accession No.	Uniprot ID:O60381
Uniprot	O60381
GeneID	26959;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HBP1

Note: This product is for in vitro research use only