

Intercellular adhesion molecule 2

Catalog No: #AP83925



Package Size: #AP83925-1 50ug #AP83925-2 100ug #AP83925-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Intercellular adhesion molecule 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ICAM-2,Lymphocyte function-associated AG-1 counter-receptor,CD102,Icam-2
Accession No.	Uniprot ID:P35330
Uniprot	P35330
GeneID	15896;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Icam2

Note: This product is for in vitro research use only