

## Proheparin-binding EGF-like growth factor

Catalog No: #AP83927



Package Size: #AP83927-1 50ug #AP83927-2 100ug #AP83927-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Proheparin-binding EGF-like growth factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	DTR,DTS,HEGFL
Accession No.	Uniprot ID:Q99075
Uniprot	Q99075
GeneID	1839;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: HBEGF

Note: This product is for in vitro research use only