

Inhibitor of nuclear factor kappa-B kinase subunit epsilon

Catalog No: #AP83928

Package Size: #AP83928-1 50ug #AP83928-2 100ug #AP83928-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inhibitor of nuclear factor kappa-B kinase subunit epsilon
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	I-kappa-B kinase epsilon,Inducible I kappa-B kinase,IKK-i,Ikke,Ikki
Accession No.	Uniprot ID:Q9R0T8
Uniprot	Q9R0T8
GeneID	56489;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ikbke

Note: This product is for in vitro research use only