

## Interleukin-2 receptor subunit beta

Catalog No: #AP83929



Package Size: #AP83929-1 50ug #AP83929-2 100ug #AP83929-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-2 receptor subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-2 receptor subunit beta,High affinity IL-2 receptor subunit beta,p70-75,CD122
Accession No.	Uniprot ID:P16297
Uniprot	P16297
GeneID	16185;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL2rb

Note: This product is for in vitro research use only