

BCLAF1 (Phospho-Ser531) Antibody

Catalog No: #12577

Package Size: #12577-1 50ul #12577-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	BCLAF1 (Phospho-Ser531) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	BCLAF1 (Phospho-Ser531) Antibody detects endogenous levels of BCLAF1 only when phosphorylated at Ser531
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human BCLAF1 (Phospho-Ser531)
Target Name	BCLAF1
Modification	Phospho
Other Names	BCLAF1, BK211L9.1, BTF, KIAA0164
Accession No.	Swiss-Prot#: Q9NYF8NCBI Gene ID: 9774
Uniprot	Q9NYF8
GeneID	9774;
Target Species	human
Calculated MW	106kd
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images

Western Blot analysis of lysates of MCF7 and MCF7 treated with LPS, using primary antibody at 1:1000 dilution.



Note: This product is for in vitro research use only