

Indolethylamine N-methyltransferase

Catalog No: #AP83931



Package Size: #AP83931-1 50ug #AP83931-2 100ug #AP83931-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Indolethylamine N-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Indolamine N-methyltransferase,Aromatic alkylamine N-methyltransferase,Amine N-methyltransferase,Thioether S-methyltransferase,TEMT,Temt
Accession No.	Uniprot ID:P40936
Uniprot	P40936
GeneID	21743;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:lnmt

Note: This product is for in vitro research use only