

Non-histone chromosomal protein HMG-17

Catalog No: #AP83935



Package Size: #AP83935-1 50ug #AP83935-2 100ug #AP83935-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Non-histone chromosomal protein HMG-17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	High mobility group nucleosome-binding domain-containing protein 2,Hmg-17,Hmg17
Accession No.	Uniprot ID:P09602
Uniprot	P09602
GeneID	15331;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmgn2

Note: This product is for in vitro research use only