

# Non-histone chromosomal protein HMG-14

Catalog No: #AP83950



Package Size: #AP83950-1 50ug #AP83950-2 100ug #AP83950-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Non-histone chromosomal protein HMG-14
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	High mobility group nucleosome-binding domain-containing protein 1,HMG14
Accession No.	Uniprot ID:P02316
Uniprot	P02316
GeneID	614915;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HMG14

Note: This product is for in vitro research use only