

Cytosolic beta-glucosidase

Catalog No: #AP83954



Package Size: #AP83954-1 50ug #AP83954-2 100ug #AP83954-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytosolic beta-glucosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytosolic beta-glucosidase-like protein 1,CBG,CBGL1
Accession No.	Uniprot ID:Q9H227
Uniprot	Q9H227
GeneID	57733;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GBA3

Note: This product is for in vitro research use only