

Growth/differentiation factor 11

Catalog No: #AP83956



Package Size: #AP83956-1 50ug #AP83956-2 100ug #AP83956-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Growth/differentiation factor 11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GDF-11, Bone morphogenetic protein 11, BMP-11, Bmp11
Accession No.	Uniprot ID:Q9Z1W4
Uniprot	Q9Z1W4
GeneID	14561;
SDS-PAGE MW	11.99
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:297-405aa Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name:Gdf11

Note: This product is for in vitro research use only