

# G-protein coupled receptor family C group 6 member A

Catalog No: #AP83958

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP83958-1 50ug #AP83958-2 100ug #AP83958-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	G-protein coupled receptor family C group 6 member A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q70VB1
Uniprot	Q70VB1
GeneID	294394;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Gprc6a

Note: This product is for in vitro research use only