

Heat shock protein HSP 90-alpha

Catalog No: #AP83960



Package Size: #AP83960-1 50ug #AP83960-2 100ug #AP83960-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock protein HSP 90-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	HSP90A,HSPCA
Accession No.	Uniprot ID:O02705
Uniprot	O02705
GeneID	397028;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSP90AA1

Note: This product is for in vitro research use only