

Histidine-rich glycoprotein

Catalog No: #AP83963



Package Size: #AP83963-1 50ug #AP83963-2 100ug #AP83963-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histidine-rich glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Histidine-proline-rich glycoprotein,HPRG,Histidine-rich glycoprotein 1,HRG1,Hrg1
Accession No.	Uniprot ID:Q99PS8
Uniprot	Q99PS8
GeneID	171016;681544;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hrg

Note: This product is for in vitro research use only