

NF-kappa-B inhibitor alpha

Catalog No: #AP83966



Package Size: #AP83966-1 50ug #AP83966-2 100ug #AP83966-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NF-kappa-B inhibitor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	I-kappa-B-alpha,IkB-alpha,Major histocompatibility complex enhancer-binding protein MAD3,IKBA,MAD3,NFKBI
Accession No.	Uniprot ID:P25963
Uniprot	P25963
GeneID	4792;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NFKBIA

Note: This product is for in vitro research use only