

## Interleukin-12 subunit beta

Catalog No: #AP83967



Package Size: #AP83967-1 50ug #AP83967-2 100ug #AP83967-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-12 subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-12B,Cytotoxic lymphocyte maturation factor 40 kDa subunit,CLMF p40,IL-12 subunit p40,NK cell stimulatory factor chain 2,NKSF2,NKSF2
Accession No.	Uniprot ID:P29460
Uniprot	P29460
GeneID	3593;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL12B

Note: This product is for in vitro research use only