

High mobility group protein B1

Catalog No: #AP83969



Package Size: #AP83969-1 50ug #AP83969-2 100ug #AP83969-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	High mobility group protein B1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	High mobility group protein 1,HMG-1,Hmg-1,Hmg1
Accession No.	Uniprot ID:P63158
Uniprot	P63158
GeneID	100862258;15289;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmgb1

Note: This product is for in vitro research use only