

Hyaluronan mediated motility receptor

Catalog No: #AP83973



Package Size: #AP83973-1 50ug #AP83973-2 100ug #AP83973-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Hyaluronan mediated motility receptor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Intracellular hyaluronic acid-binding protein, Receptor for hyaluronan-mediated motility, CD168, Ihabp, Rhamm |
| Accession No. | Uniprot ID:Q00547 |
| Uniprot | Q00547 |
| GeneID | 15366; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hmnr

Note: This product is for in vitro research use only