

Interferon epsilon

Catalog No: #AP83974



Package Size: #AP83974-1 50ug #AP83974-2 100ug #AP83974-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon epsilon
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	IFN-epsilon,Interferon epsilon-1,IFNE1
Accession No.	Uniprot ID:A7UHZ5
Uniprot	A7UHZ5
GeneID	100125969;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IFNE

Note: This product is for in vitro research use only