

Heat shock protein beta-6

Catalog No: #AP83980



Package Size: #AP83980-1 50ug #AP83980-2 100ug #AP83980-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Heat shock protein beta-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	HspB6
Accession No.	Uniprot ID:Q148F8
Uniprot	Q148F8
GeneID	534551;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HSPB6

Note: This product is for in vitro research use only