

## Serine protease HTRA2, mitochondrial

Catalog No: #AP83981



Package Size: #AP83981-1 50ug #AP83981-2 100ug #AP83981-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine protease HTRA2, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	High temperature requirement protein A2,HtrA2,Omi stress-regulated endoprotease,Serine protease 25,Serine proteinase OMI,OMI,PRSS25
Accession No.	Uniprot ID:O43464
Uniprot	O43464
GeneID	27429;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:HTRA2

Note: This product is for in vitro research use only