

Hematopoietic SH2 domain-containing protein

Catalog No: #AP83986



Package Size: #AP83986-1 50ug #AP83986-2 100ug #AP83986-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hematopoietic SH2 domain-containing protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Hematopoietic SH2 protein, Adaptor in lymphocytes of unknown function X, Alx
Accession No.	Uniprot ID:Q6VYH9
Uniprot	Q6VYH9
GeneID	209488;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hsh2d

Note: This product is for in vitro research use only