

Hyaluronidase-1

Catalog No: #AP83988



Package Size: #AP83988-1 50ug #AP83988-2 100ug #AP83988-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hyaluronidase-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Hyal-1,Hyaluronoglucosaminidase-1
Accession No.	Uniprot ID:Q6RHW4
Uniprot	Q6RHW4
GeneID	404698;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HYAL1

Note: This product is for in vitro research use only