

Heat shock cognate 71 kDa protein

Catalog No: #AP83989



Package Size: #AP83989-1 50ug #AP83989-2 100ug #AP83989-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock cognate 71 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Heat shock 70 kDa protein 8,Hsc70,Hsc73
Accession No.	Uniprot ID:P63018
Uniprot	P63018
GeneID	24468;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspa8

Note: This product is for in vitro research use only