

E3 ubiquitin-protein ligase HUWE1

Catalog No: #AP83990



Package Size: #AP83990-1 50ug #AP83990-2 100ug #AP83990-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase HUWE1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HECT, UBA and WWE domain-containing protein 1,HECT-type E3 ubiquitin transferase HUWE1,Upstream regulatory element-binding protein 1,URE-B1,Ureb1
Accession No.	Uniprot ID:P51593
Uniprot	P51593
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Huwe1

Note: This product is for in vitro research use only