

Keratin, type I cytoskeletal 10

Catalog No: #AP83992



Package Size: #AP83992-1 50ug #AP83992-2 100ug #AP83992-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokeratin-10,CK-10,Keratin-10,K10,Type I keratin Ka10,Ka10
Accession No.	Uniprot ID:Q6IFW6
Uniprot	Q6IFW6
GeneID	450225;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Krt10

Note: This product is for in vitro research use only