

Keratin, type I cytoskeletal 18

Catalog No: #AP83993



Package Size: #AP83993-1 50ug #AP83993-2 100ug #AP83993-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokeratin-18,CK-18,Keratin-18,K18,Krt1-18
Accession No.	Uniprot ID:Q5BJY9
Uniprot	Q5BJY9
GeneID	294853;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Krt18

Note: This product is for in vitro research use only