

Hemopexin

Catalog No: #AP83994



Package Size: #AP83994-1 50ug #AP83994-2 100ug #AP83994-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Hemopexin |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Pig |
| Other Names | Hyaluronidase |
| Accession No. | Uniprot ID:P50828 |
| Uniprot | P50828 |
| GeneID | 396998; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HPX

Note: This product is for in vitro research use only