

Hydroxymethylglutaryl-CoA synthase, cytoplasmic

Catalog No: #AP83995



Package Size: #AP83995-1 50ug #AP83995-2 100ug #AP83995-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hydroxymethylglutaryl-CoA synthase, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HMG-CoA synthase,3-hydroxy-3-methylglutaryl coenzyme A synthase
Accession No.	Uniprot ID:Q8JZK9
Uniprot	Q8JZK9
GeneID	208715;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmgcs1

Note: This product is for in vitro research use only