

# 3-hydroxy-3-methylglutaryl-coenzyme A reductase

Catalog No: #AP83996



Package Size: #AP83996-1 50ug #AP83996-2 100ug #AP83996-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	3-hydroxy-3-methylglutaryl-coenzyme A reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	HMG-CoA reductase, HMGR
Accession No.	Uniprot ID:Q1W675
Uniprot	Q1W675
GeneID	100144446;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: HMGR

Note: This product is for in vitro research use only