

## Interleukin-17 receptor B

Catalog No: #AP84006



Package Size: #AP84006-1 50ug #AP84006-2 100ug #AP84006-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-17 receptor B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-17 receptor B,Cytokine receptor-like 4,IL-17 receptor homolog 1,IL-17Rh1,Interleukin-17B receptor,IL-17B receptor,CRL4,EVI27,IL17BR,UNQ2501/PRO19612
Accession No.	Uniprot ID:Q9NRM6
Uniprot	Q9NRM6
GeneID	55540;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL17RB

Note: This product is for in vitro research use only