

Insulin-like growth factor-binding protein 3

Catalog No: #AP84008



Package Size: #AP84008-1 50ug #AP84008-2 100ug #AP84008-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor-binding protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IBP-3
Accession No.	Uniprot ID:P20959
Uniprot	P20959
GeneID	282261;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IGFBP3

Note: This product is for in vitro research use only