

## Interleukin-17B

Catalog No: #AP84010



Package Size: #AP84010-1 50ug #AP84010-2 100ug #AP84010-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-17B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-17B,Cytokine Zcyto7,Interleukin-20,IL-20,Neuronal interleukin-17-related factor,IL20,NIRF,ZCYTO7,UNQ516/PRO1031
Accession No.	Uniprot ID:Q9UHF5
Uniprot	Q9UHF5
GeneID	27190;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL17B

Note: This product is for in vitro research use only