

# Interferon gamma receptor 1

Catalog No: #AP84012



Package Size: #AP84012-1 50ug #AP84012-2 100ug #AP84012-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interferon gamma receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IFN-gamma receptor 1, Interferon gamma receptor alpha-chain, IFN-gamma-R-alpha, CD119, Ifngr
Accession No.	Uniprot ID:P15261
Uniprot	P15261
GeneID	15979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Ifngr1

Note: This product is for in vitro research use only