

Heat shock factor protein 4

Catalog No: #AP84014

Package Size: #AP84014-1 50ug #AP84014-2 100ug #AP84014-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock factor protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HSF 4,Heat shock transcription factor 4,HSTF 4
Accession No.	Uniprot ID:Q9ULV5
Uniprot	Q9ULV5
GeneID	3299;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSF4

Note: This product is for in vitro research use only