

Furin

Catalog No: #AP84015



Package Size: #AP84015-1 50ug #AP84015-2 100ug #AP84015-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Furin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Dibasic-processing enzyme, Paired basic amino acid residue-cleaving enzyme, PACE, Prohormone convertase 3, Fur, Pcsk3
Accession No.	Uniprot ID:P23188
Uniprot	P23188
GeneID	18550;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Furin

Note: This product is for in vitro research use only