

Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-3

Catalog No: #AP84018

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84018-1 50ug #AP84018-2 100ug #AP84018-3 1mg

Description

Product Name	Guanine nucleotide-binding protein G(I)/G(S)/G(T) subunit beta-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Transducin beta chain 3
Accession No.	Uniprot ID:P16520
Uniprot	P16520
GeneID	2784;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: GNB3

Note: This product is for in vitro research use only