

Interleukin-4

Catalog No: #AP84020



Package Size: #AP84020-1 50ug #AP84020-2 100ug #AP84020-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-4,B-cell IgG differentiation factor,B-cell growth factor 1,B-cell stimulatory factor 1,BSF-1,IGG1 induction factor,Lymphocyte stimulatory factor 1,IL-4
Accession No.	Uniprot ID:P07750
Uniprot	P07750
GeneID	16189;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL4

Note: This product is for in vitro research use only