

# Solute carrier family 2, facilitated glucose transporter member 3

Catalog No: #AP84025

Package Size: #AP84025-1 50ug #AP84025-2 100ug #AP84025-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Solute carrier family 2, facilitated glucose transporter member 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	CEF-GT3,Glucose transporter type 3,GLUT-3,GLUT3
Accession No.	Uniprot ID:P28568
Uniprot	P28568
GeneID	396517;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC2A3

Note: This product is for in vitro research use only