

## Hemoglobin subunit beta-1

Catalog No: #AP84027



Package Size: #AP84027-1 50ug #AP84027-2 100ug #AP84027-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hemoglobin subunit beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Beta-1-globin,Hemoglobin beta-1 chain,Hemoglobin beta-major chain
Accession No.	Uniprot ID:P02088
Uniprot	P02088
GeneID	100503605;101488143;15129;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Hbb-b1

Note: This product is for in vitro research use only