

High mobility group protein B3

Catalog No: #AP84028



Package Size: #AP84028-1 50ug #AP84028-2 100ug #AP84028-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	High mobility group protein B3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q32L31
Uniprot	Q32L31
GeneID	613729;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HMGB3

Note: This product is for in vitro research use only