

Isocitrate dehydrogenase [NADP] cytoplasmic

Catalog No: #AP84033



Package Size: #AP84033-1 50ug #AP84033-2 100ug #AP84033-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Isocitrate dehydrogenase [NADP] cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IDH,Cytosolic NADP-isocitrate dehydrogenase,IDP,NADP(+)-specific ICDH,Oxalosuccinate decarboxylase,ICDH
Accession No.	Uniprot ID:Q9XSG3
Uniprot	Q9XSG3
GeneID	281235;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IDH1

Note: This product is for in vitro research use only