

## Interferon alpha-inducible protein 6

Catalog No: #AP84034



Package Size: #AP84034-1 50ug #AP84034-2 100ug #AP84034-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interferon alpha-inducible protein 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Interferon-induced protein 6-16, Ifi-6-16
Accession No.	Uniprot ID:Q6IED8
Uniprot	Q6IED8
GeneID	512913;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: IFI6

Note: This product is for in vitro research use only