

Heat shock protein beta-2

Catalog No: #AP84036

Package Size: #AP84036-1 50ug #AP84036-2 100ug #AP84036-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock protein beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	HspB2,DMPK-binding protein,MKBP
Accession No.	Uniprot ID:Q16082
Uniprot	Q16082
GeneID	3316;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HSPB2

Note: This product is for in vitro research use only