

Histone acetyltransferase type B catalytic subunit

Catalog No: #AP84037



Package Size: #AP84037-1 50ug #AP84037-2 100ug #AP84037-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone acetyltransferase type B catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Histone acetyltransferase 1,KAT1
Accession No.	Uniprot ID:O14929
Uniprot	O14929
GeneID	8520;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:HAT1

Note: This product is for in vitro research use only